



Escape-Mattress Bariatric



For the evacuation of obese individuals

What's an Escape Mattress?

The **Escape-Mattress Bariatric** is designed for evacuating obese people (up to 300kg) in lying position. The **Escape-Mattress Bariatric** needs to be operated by several trained persons. Also suitable for spiral staircases.

The **Escape-Mattress Bariatric** can be used to move evacuees along the floor (by using at least 4 persons and using the additional towing straps) and have 12 handles placed on the sides to allow lifting and easy maneuverability.

Advantages of the EscapeMattress Compact:

- TÜV certified on 300KG.
- Extra wide version.
- Made of high quality materials (fire resistance classification M2).
- Extra-long towing straps to help operators with improved handling practices.
- Handgrips to assist with sideways manoeuvrability.
- Carrying people is no longer necessary.
- Ready for use within seconds.
- Safe evacuation down the stairs.
- Speed of descent is determined by the user.
- Fits into your emergency evacuation plan.
- CE certified.

Specifications Escape Mattress Bariatric

Article number:	21080
Dimensions (length x width x thick):	± 200 x 86 x 2 cm
Maximum payload:	300 kg
Weight:	±10 kg incl. bag
Demonstration and advice provided:	✓
Equipped with additional, movable drag straps:	✓
Ready for use within seconds:	✓
Speed of descent to be determined by the user:	✓
Blood-, urine-, acid- & alkali resistant:	✓
Fire resistance classification M2:	✓
Equipped with headrest, body-covers & foot sack:	✓
Including storage bag:	✓
Warranty:	2 years

Our advice:

Make sure you have an **Escape-Mattress Bariatric** available, so you can help obese people who rely on you to be rescued safely. Even when you need to tackle spiral staircases.

"Failure to prepare for an emergency has its consequences."



Escape Mobility Middle East.

NAFFCO FFZ

World Headquarters, Jebel Ali Free Zone, Dubai, U.A.E., P.O. Box: 26216

W www.escape-mobility.ae E info@naffco.com

T +971 4 815 1111

F +971 4 815 1222

