



Escape-Sheet Standard



Best choice for hospitals and healthcare facilities

What's an Escape Sheet?

People who are unable to leave their beds, must be evacuated during an emergency. One way of evacuation is an **Escape-Sheet**. An **Escape-Sheet** is permanently under the mattress and in an event of an emergency, the safety straps of the **Escape-Sheet** are secured over the mattress and the patient. The patient is then secure and immobilized for safe horizontal and vertical evacuation down the stairs.

The **Escape-Sheet** is an "A-label quality" product and produced in-house.

Our sheets are CE certified & TÜV certified up to 250kg. We guarantee the highest quality standards.

The **Escape-Sheet** is positioned permanently between the mattress and the bed frame. In the event of an emergency, the safety straps of the **Escape-Sheet** are secured over the mattress and the patient.

The patient is immobilised safely for transporting down the stairs. The **Escape-Sheet** is used frequently in hospitals and nursing homes.

The **Escape-Sheet** is widely used by hospitals and healthcare facilities. Assists in moving a patient to a place of safety, horizontally as well as vertically. Our sheets are water repellent, anti-static and also available in heat and flame resistant.

Advantages of the Escape Sheet Standard:

- Straps and handles highly visible in bright green colour.
- Extra-long towing straps.
- Plastic buckle pockets increase the life expectancy of the product
- Handgrips to assist with sideways manoeuvrability.
- Static-Control™ (anti-static polyester).
- Adjusting corner straps (suitable for each mattress height).
- Hydro-tech Nano technology (water-repellent)

Specifications Escape Sheet Standard	
Article number:	21240
Dimensions (length x width):	200 x 90 cm
Material:	Polyester
Maximum payload:	250 kg
Weight:	1000 gram
Demonstration and advise provided:	✓
Ready to use within a few seconds:	✓
Immediate evacuation directly from the bed:	✓
Industrial wash 75° C:	✓
Hydro-tech nano technology (water-repellent):	✓
Adjustable corner straps:	✓
For all mattress thicknesses:	< 23 cm
Buckles in pockets:	✓
Heat and flame resistant version (HFR):	Optional
Towing straps available in 3 different colors:	Optional
Warranty:	2 Years

Usage:

The **Escape-Sheet** is widely used by hospitals and healthcare facilities. Assists in moving a patient to a place of safety.



Escape Mobility Middle East.
 NAFFCO FFZ
 World Headquarters, Jebel Ali Free Zone, Dubai, U.A.E., P.O. Box: 26216
W www.escape-mobility.ae **E** info@naffco.com
T +971 4 815 1111 **F** +971 4 815 1222

